



ProCair Prime

The ground-breaking, unique alternating pressure mattress with an integrated pump.



Unique design reduces facilities' repair cost by integrating pump into mattress



Can be adapted between dynamic mode or hybrid mode based on different clinical needs



Cord retention management to prevent tripping hazard

ProCair Prime

Easy to setup, the ProCair Prime automatically adjusts the pressure to suit the weight of the patient and does not need a separate pump unit to hang from the bed - it is all built in!



DYNAMIC MODE

HYBRID MODE



The ProCair Prime is adaptable for different clinical needs. It can be used as a fully alternating pressure mattress (dynamic mode) or easily converted to hybrid mode by positioning the foam sheet above the air cell layer.

SPECIFICATIONS

CODE	APMPC-R03	APMPC-R03K
CAPACITY (SAFE WORKING LOAD)	30–200 kg	30–200 kg
SYSTEM ALTERNATION	1 in 2	1 in 2
WARRANTY	3 year	3 year
PRESSURE INJURY RISK	Very high risk	Very high risk
FOWLER BOOST	N/A	Manual
ARTG	289458	289458
MATTRESS DIMENSIONS (L x W x H)	2000 x 900 x 180mm	2000 x 1050 x 180mm
MATTRESS WEIGHT	12 kg	13 kg



Refer to the Technical Data Sheet for detailed specifications.

