Shipping adress ☐ see adress

☐ Order
☐ Estimate

Purchaser: _____

Telephone: _____

MADE IN GERMANY
CRASH TESTED ISO 7176-19
UNITOMIC® ANATOMIC SEAT&BACK-FRAME
MATRIX ANATOMIC SEAT&BACK-SYSTEM

Configuration: **BASIC-Line**

Little

Kids & Teens

Activ children wheelchair made of aluminium, foldable
Frame straight or with abduction

All dimensions in cm
Dimensions tol. +/- 0,5cm

Capacity up to
100 kg (SW 24 and 26cm)
120 kg (SW 28 till 36cm)

user weight: _____ kg


Little SLT frame Little AR

Frame with **integrated** footrest (without footrest and without side-guard)

frame straight



SD 26 - 38



Abduction frame *




SD 26 - 32 frame short
SD 34 - 38 frame long

* from 4° camber possible

☐ Code 2421 ☐ Code 2422

Length of the frame



frame straight
51cm = standard (SD 26-38)

Abduction frame
54cm = short (SD 26-32)
57cm = long (SD 34-38)


- Total length depends on the wheel position and position of the footrest
- Frame length for UNITOMIC L+3cm

Frame with **detachable footrest** (without footrest and without side-guard)

frame straight



SD 26 - 34 frame middle
SD 36 - 38 frame long



Abduction frame *




SD 26 - 34 frame short
SD 36 - 38 frame long

* from 4° camber possible

☐ Code 2428 ☐ Code 2429

Length of the frame with foot-rest holder



frame straight
54cm = middle (SD 26-34)
58cm = standard (SD 36/38)

Abduction frame
52cm = short (SD 26-34)
56cm = long (SD 36/38)

Technical changes and print mistakes are reserved.

basic measures

Seat width (SW)

☐ 24 ☐ 26 ☐ 28 ☐ 30
☐ 32 ☐ 34 ☐ 36

Seat depth (SD)

☐ 26 ☐ 28 ☐ 30 ☐ 32
☐ 34 ☐ 36 ☐ 38

Please note the details for frame long and short!!

Back height (BH)

☐ 25 ☐ 30
☐ 35 ☐ 40

(+2,5cm adjustable)

Lower Leg Length (LLL):

LLL: ..... cm
(regarding seat upholstery)

min. 12 cm
max. 45 cm (without cushion)

max. LLL = SH - 3cm

Technical data see back side

* frame straight Code 2421 not possible (only with abduction frame Code 2422) !!!

seat height





Seat height in cm
(without cushion)

front: _____ cm
(Standard = 48/49cm)

back: _____ cm
(Standard = 45/46cm)


VD2362E-29

DOAbility

Technical data


seat width **SW**

from outside crossbar to outside crossbar +6mm




seat depth **SD**

from back tube to upper edge seat upholstery



back height **BH**

from upper edge seat upholstery to upper edge back upholstery



Lower Leg Length LLL:		15 to 42 cm
Total width wheelchair:	with wheels and handrim	SW + 18 cm
	with drum brake	SW + 20 cm
Folding measurement:	with wheels	about 30cm
	with drum brake	about 32cm
	without wheels	about 22cm

Seat height **SH**

Little

Rear wheel	Castor fork	minimum seat height front	minimum seat height back	maximum seat height front	maximum seat height back
20 inch	4 inch - short	44 cm	36 cm*	46 cm	46 cm
	5 inch - short	44 cm	36 cm*	47 cm	46 cm
	6 inch - short	45 cm	36 cm*	48 cm	46 cm
	140 mm - wide	46 cm	36 cm*	50 cm	46 cm
22 inch	4 inch - short	43 cm	39 cm*	46 cm	46 cm
	4 inch - long	46 cm	39 cm*	50 cm	49cm
	5 inch - short	44 cm	39 cm*	47 cm	46 cm
	5 inch - long	47 cm	39 cm*	51 cm	49 cm
	6 inch - short	45 cm	39 cm*	48 cm	48 cm
	6 inch - long	48 cm	39 cm*	52 cm	51 cm
	140 mm - wide	46 cm	39 cm*	50 cm	49 cm

All measurements ± 0,5 cm

*Advice: with drum brake you have to add 2cm!!

Notice: at the castor wheels PU-black are the minimal seat height with the fork short and extra short not possible

Technical changes and print mistakes are reserved.

DOAbility

frame colour

Because of the material are different colours possible

shiny

- ☐ 25 RAL9010 white
- ☐ 56 RAL1021 goldyellow met.
- ☐ 24 RAL1028 melonyellow
- ☐ 27 mandarin met.
- ☐ 41 RAL3020 traffic red met.
- ☐ 06 dormantred
- ☐ 40 redmetallic
- ☐ 43 RAL4010 telemagenta met.
- ☐ 53 cassis
- ☐ 59 atlanticblue met.
- ☐ 44 RAL5015 sky blue met.
- ☐ 16 RAL5002 ultramarine blue
- ☐ 17 dormantblue
- ☐ 64 petrolgreen met.
- ☐ 50 applegreen
- ☐ 51 applegreen metallic
- ☐ 45 RAL6018 yellowgreen
- ☐ 60 mojitogreen met.
- ☐ 19 silvermetallic
- ☐ 23 RAL9005 black

met. = metalized

frame colour

Because of the material are different colours possible

mat / satined

- ☐ 76 white satined
- ☐ 65 goldyellow met. satined
- ☐ 66 mandarin met. satined
- ☐ 67 traffic red satined
- ☐ 68 telemagenta met. satined
- ☐ 69 purpur met. satined
- ☐ 61 cool mint mat
- ☐ 72 skyblue met. satined
- ☐ 71 atlanticblue met. satined
- ☐ 70 ultramarine blue satined
- ☐ 75 petrolgreen met. satined
- ☐ 73 yellowgreen met. satined
- ☐ 74 mojitogreen met. satined
- ☐ 77 silvermetallic satined
- ☐ 22 RAL9011 black mat

met. = metalized

upholstery colours

upholstery standard

- ☐ DE290 black





upholstery breathable fabric

- ☐ AG410 yellow
- ☐ AG430 red
- ☐ AG441 pink
- ☐ AG450 ultramarine
- ☐ AG451 king blue
- ☐ AG460 lemon
- ☐ AG470 grey
- ☐ AG490 black

only for adjustable back upholstery

Code 201

frame accessories

Tip shackle	Transit wheels	Walking stick holder, aluminium, with a loop	Anti-tip wheel	Anti-tip wheel high adjustabe and with gear segment	
☐ right ☐ left 		for 2 crutches ☐ right ☐ left 	☐ right ☐ left 	☐ right ☐ left 	
Code 105	Code 102	Code 103	Code 104	Code 2405	
☐ Piece	☐ Pair	☐ Piece	☐ Piece	☐ Piece	

note

Technical changes and print mistakes are reserved.

DOAbility

Edition-Line Little	BASIC-Line	TREND-Line	HIGH-Line
	The **BASIC-Line** is the standard equipment. There are only certain parts in the frame colours powder-coated and the upholstery is standard. ☐ Code 439B	If you order the **TREND-Line** there are a variety of attachments also in colour of the wheelchair. And it is possible, to select a breathable upholstery. **Code 439T**	The **HIGH-Line** optimizes the design by ordering many attachments in colour of the wheelchair. And it is possible to order the upholstery in SLIMLINE-Design. **Code 439H**
	standard	☐ Piece	☐ Piece


frame accessories

Overview for **Edition-Line** Little

	BASIC-Line	TREND-Line	HIGH-Line
axle plate	black anodized		
shell for axle plate		black anodized	black anodized
handrim aluminium			
castor adapter			
cap for castor adapter			colour of the wheelchair
cap for extender		colour of the wheelchair	
fork			
footplate	colour of the wheelchair		
front frame parts	colour of the wheelchair	colour of the wheelchair	colour of the wheelchair
back frame parts			
cross bar			
footrest holder (Little 2428/2429)			
side guard			
upholstery standard	upholstery	upholstery	upholstery
upholstery breathable fabric and uni coloured			
upholstery breathable fabric SLIMLINE-Design			
Edition-Line Little	Standard	include CODE 201	include CODE 201 include CODE 202 **SlimLine**-Design for adjustable, breathable back upholstery 1 2 bicoulored upholstery

(other data see Page 3)



frame

Frame part +30mm backwards only for straight back	Spasticity frame back frame by street (reinforced version)	
30mm Code 204  back tube 30mm displaced backwards (e.g. for straight back with back shell)	recommended with strong spasticity Code 426	
☐	☐	





Technical changes and print mistakes are reserved.

DOAbility




footrest

Single footrest, angle adjustable from LLL 35	Single footrest, angle adjustable and mounted in frame LLL 15 - 32cm	
 Code 2410	 Code 2411*	* Note: Frame straight from SW 28 Abduction frame from SW 24
☐ Piece	☐ Piece	


footrest

Single footrest **alu-profile**, angle adjustable, in colour from LLL 35cm	Single footrest **flat-profile**, angle adjustable, in colour from LLL 35cm	Single footrest **alu-profile**, mounted in frame with clamp plate from LLL 12 - 35 from SW 26	Single footrest **flat-profile**, mounted in frame with clamp plate from LLL 12 - 35 from SW 26	
 Code 2412	 Code 2432	 Clamp plate: ☐ 60mm ☐ 85mm ☐ 145mm Code 2413	 Clamp plate: ☐ 60mm ☐ 85mm ☐ 145mm Code 2433	
☐ Piece	☐ Piece	☐ Piece	☐ Piece	

footrest

Double footrest **alu-profile**, angle adjustable, in colour from LLL 35 from SW 26	Double footrest **flat-profile**, angle adjustable, in colour from LLL 35 from SW 26	Double footrest **alu-profile**, mounted in frame with clamp plate from LLL 12 - 35 from SW 28	Double footrest **flat-profile**, mounted in frame with clamp plate from LLL 12 - 35 from SW 28	
 Code 2414	 Code 2434	 Clamp plate: ☐ 60mm ☐ 85mm ☐ 145mm Code 2415	 Clamp plate: ☐ 60mm ☐ 85mm ☐ 145mm Code 2435	
☐ Pair	☐ Pair	☐ Piece	☐ Piece	

footrest

Single footrest **alu-profile, displaced** with footrest tube and clamp plate, angle adjustable from LLL 28	Double footrest **alu-profile, displaced** with footrest tube and clamp plate, angle adjustable from LLL 28 from SW 28		Spasticity locking for **single** footrest footrest **alu-profile** Code 122 ☐ | footrest **alu. flat-profile** Code 1122 ☐	Spasticity locking for **double** footrest footrest **alu-profile** Code 123 ☐ | footrest **alu. flat-profile** Code 1123 ☐
 Clamp plate 85mm Code 2417	 Clamp plate 60/85mm Code 2416			
☐ Piece	☐ Piece		☐ Piece	☐ Pair

note





Technical changes and print mistakes are reserved.

DOAbility

Height of push handles

Height of push handles for companion = Back height + Seat height + Height of push handle

Height of push handles for companion = BH + SH + HPH




	Back standard	Back with height adjustable push handle				
		15cm normal	35cm long	45cm extra long reinforced	cropped	bent rewards
						
HPH						
min.	3cm	4cm	4cm	4cm	21cm	12cm
max.		19cm	39cm	49cm	38cm	29cm
					push handle each side 6cm outside	push handle each side 9cm backwards

Technical changes and print mistakes are reserved.

Berollka-aktiv Rollstuhltechnik GmbH, [illegible] 74889 Sinsheim, Phone +49 7261-7351-0, Fax +49 7261-7351-10
E-Mail: info@berollka.de [illegible]

DOAbility
Our printed terms and conditions apply to all orders.

back

Back with adjustable upholstery	Back with adjustable upholstery and hight adjustable handles	Gear segment for back angle adjustion for **straight** back		SlimLine-Design for adjustable, breathable back upholstery (fabric 74, 75, 89, 92, 94, 95 und 96)
	☐ 15cm (standard) ☐ 35cm (long) ☐ cropped ☐ bent rewards 	from BH 35  Note: Only in combinaton with CODE 447 or 448 possible !!		 1 2 Colour 1 (inside): _____ Colour 2 (outside): _____
☐ Code 447	Code 448	Code 407		bicoloured upholstery
Standard	☐ Piece	☐ Piece		☐ Code 202

back

* Attention: with UNITOMIC back choose seat depth 2cm shorter

Back **UNITOMIC** * with adjustable upholstery (anatomical back)	Back **UNITOMIC** * with adjustable upholstry and adjustable handles	Gear segment for back angle adjustion for **UNITOMIC** back	**FIDIBUS back with clips**
from BH 35  Note: Choose seat-depth 2cm shorter	from BH 35  Note: Choose seat-depth 2cm shorter	from BH 40  Note: only in combination with CODE 410 or 412 possible !!	ergonomic **back cushion with clips** extrem light ergonomic formed back  curved back 30mm (standard) ☐ Code 2292F curved back 50mm ☐ Code 2292F-5 Please take prices and further information from FIDIBUS order form price reduction of the back cover is consider in the price of the FIDIBUS back
Code 412	Code 410	Code 411	
☐ Piece	☐ Piece	☐ Piece	☐ Piece


Note: When using SlimLine-Design (Code 202) the seat cushion (Code 470 / 471 / 472) is also supplied in SlimLine-Design.

back

FIDIBUS ASR back

anatomical shaped back plate ASR with Clips	anatomical shaped back plate ASR lockable	anatomical shaped back cushion ASR
Alu-back plate **coloured**, curved and beveled (without cushion) colour: frame colour	Alu-back plate **coloured**, cuvred and beleved (without cushion) colour: frame colour	back cushion (foam) with cover (without plate)
CODE	CODE	CODE
30mm curved back (standard) ☐ 2294F	30mm curved back (standard) ☐ 2294FV	30mm curved back (standard) ☐ 2296A
50mm curved back ☐ 2294F-5	50mm curved back ☐ 2294FV-5	50mm curved back ☐ 2296A-5
☐ Piece	☐ Piece	☐ Piece





seat

Cushion compounded foam with black cover	Cushion, compounded foam with breathable cover	Change cover for cushion
		Standard-fabric black ☐ Code 121
Code 119 ☐ 3 cm	Code 471 ☐	**breathable** fabric ☐ Code 470
Code 120 ☐ 5 cm	Code 472 ☐	Colour: _____
Piece		Piece

seat

Seat belt black	Safty seat belt with cushion, black		Seat cushion by wheelchair with back upholstery	Seat cushion by wheelchait with MATRIX- or FIDIBUS back
Plastic closure Code 227 ☐ Piece	2-point-lap-sash seat belt with metal lock		FIDIBUS-F ergonomic seat cushion **straight at the back** extrem light ergonomic formed seat cushion from SW 24 from SD 26  Code 2291F Please take prices and further informations from the FIDIBUS oder form	FIDIBUS ergonomic seat cushion extrem light ergonomic formed seat cushion from SW 24 from SD 26  Code 2291 Please take prices and further informations from the FIDIBUS oder form
Metal closure Code 125	Code 124			
☐ Piece	☐ Piece		☐ Piece	☐ Piece





accessories

Head rest adjustable (for back with CODE 447 and 412)	Head rest adjustable (for height adjustable handles) (for back with CODE 448, and 410)	Holder for head rest for height adjustable handles for back with CODE 448 and 410		Stabilization bar only for height adjustable handles	Stabilization bar push handles (= pushing bar)
		for head rest tube: ☐ 15x15mm Code 225V ☐ Ø 19mm Code 225R 			
Code 224	Code 225			Code 178	Code 478
☐ Piece	☐ Piece	☐ Piece		☐ Piece	☐ Piece



Technical changes and print mistakes are reserved.

DOAbility




side guards

Clothing guard **round** screwed on, **aluminium** with cold protection (in frame color)  Code 433	Clothing guard screwed on, **plastic,** black  Code 430	Clothing guard **round PC-Carbon,** with cold protection, screwed on (CARBON-Optics)  Code 230	Clothing side guard **ERGO AL** aluminium with cold protection, screwed on (in frame colour)  Code 437	
☐ Pair	☐ Pair	☐ Pair	☐ Pair	



side guards

Clothing guard **round** pluggable, **aluminium** with cold protection (in frame color)  Code 434	Clothing guard pluggable, **plastic,** black  Code 431	
☐ Pair	☐ Pair	




side guards

Clothing guard **round** screwed on, **aluminium** with cold protection, armrest screwed on (side guard in frame colour) armrest short  Code 231	Clothing side guard **ERGO AL** Aluminium with armrest short screwed on, with cold protection (side guard in frame colour) armrest short  Code 4371		Arm pads swingable to side and high adjustable ☐ armrest short ☐ armrest long  Code 432 (without clothing guard)
☐ Pair	☐ Pair		☐ Pair

side guards

Side guard screwed on, aluminium, with cold protection  Code 232	Side guard with armrest short, screwed on, aluminium, armrest short armrest short  Code 238	
☐ Pair	☐ Pair	

therapy table

Therapy table clamped in the frame with clamp bracket at the frame  Code 138	Therapy table **plug** on the frame with plug attachment on the frame  Code 139	Therapy table on the frame quick release clamped (at the frame)  Code 138S	
☐ Piece	☐ Piece	☐ Piece	

Technical changes and print mistakes are reserved.

DOAbility

Combinationen wheels with frame length

rear wheel	Little SLT (frame integrated footrest) CODE 2421 frame straight (integrated)	Little SLT (frame integrated footrest) CODE 2422 frame with abuction (integrated)	Little AR (with demounted footrest) CODE 2428 frame straight (with demountable fottrest)	Little AR (with demounted footrest) CODE 2429 frame with abdukction (with demountable footrest)
20"	possible	possible	possible only with 4"/5"	possible only with 4"/5"
22"	possible	possible	possible only with 4"/5"	ST 26 - 34 Frame short possible (only with castor 4" or 5" possible) (only wheel position „C" possible!!!) ST 36 - 38 Frame long possible (only castor 4" or 5" possible) (only wheel position „B" and „C" possible!!!)

Wheel plate:

Center of gravity adjustment (measured from black tube to wheel center)


Wheel plate staggered / back straight:

A = 8,0 cm

B = 6,0 cm

C = 4,0 cm

wheel plate turned / backwards B = 2,0 cm



-> if no information from the customer, we mount the wheel plate at position C
(position B at ERGO-AL)
position A only possible with side guard ERGO-AL Code 437 !!)


Wheel plate staggered / UNITMOIC back:

A = 8,0 cm (3,5cm at bow)

B = 6,0 cm (1,5cm at bow)

C = 4,0 cm (-0,5cm at bow)

wheel plate turned / backwards B = 2,0 cm (-2,5cm at bow)



-> if no information from the customer, we mount the wheel plate at position B

Technical changes and print mistakes are reserved.

DOAbility

Note: the combination 22" and 6" at short frame (Code 2428 or 2429) is not possible !!!!

rear wheel

Rear wheel standard (alu-rim) with quick release axle

Pneumatic			Solid		
☐ 20"x1*	fine profile	Code 142	☐ 20"x1 3/8 grey		CODE 143
☐ 20"x1	black fine profile	Code 142S	☐ 20"x1 grey		CODE 240
☐ 22"x1 3/8*	rough profile	Code 150	☐ 20"x1 black		CODE 241
☐ 22"x1*	rough profile	Code 147	☐ 22"x1 3/8 grey		CODE 149
☐ 22"x1*	fine profile	Code 144	☐ 22"x1 grey		CODE 141
☐ 22"x1	black fine profile	Code 144S	☐ 22"x1 black		CODE 141S

*grey tyres

Standard

22" x1 Performance tyre

☐ 22"x1 Code 577

Marathon Plus Evolution (black pneumatic tyre)

Right Run **22"x1** Code 576 ☐ black

Camber
- 0° ☐
- 2° ☐
- 4° ☐ (standard)
- 7° ☐
- 10° ☐

handrim for rear wheel

Handrim for rear wheel

aluminium **black**	aluminium **silver**	stainless steel	TETRA-GRIP (Super Grip)	Maxgrepp ERGO Para	Maxgrepp ERGO Para CP	CURVE black anodized
☐ 20" Code 158		☐ 20" Code 291	☐ 20" Code 254	☐ 20" Code 558	☐ 20" Code 559	
☐ 22" Code 455	☐ 22" Code 492	☐ 22" Code 155	☐ 22" Code 252	☐ 22" Code 550	☐ 22" Code 551	☐ 22" Code 4502
extension narrow ☐ wide ☐			extension narrow	extension narrow	extension narrow	extension narrow

Standard

Handrim cover black

☐ 20 inch Code 290

☐ 22 inch Code 191

☐ Pair

rear wheel

Rear wheel with **drum brake** and quick release axle (for companion)

☐ 20 inch Code 2487

☐ 22 inch Code 2489

☐ Pair

light weight wheel

Light weight wheel **LWW** with black rim, 24 radial spokes black with alu hand rim and pneumatic tyre (fine profile)


	grey tyre		black tyre
	hand rim **black**	hand rim **silver**	hand rim **black**
20 x1"	☐ CODE 565	☐ CODE 565S	☐ CODE 2868
22 x1"	☐ CODE 564	☐ CODE 564S	☐ CODE 2867

☐ Pair not possible with solid tyre !

wheel plate for rear wheel


Center of gravity adjustment

wheel position: _____



staggered wheel plate

Wheel extension for drum brake




Code 162

only with wheel plate straight

☐ Pair

Wheel extension



Code 161

only with wheel plate straight

☐ Pair

rear wheel - spoke cover

Spoke cover transparent



☐ 20 inch Code 164

☐ 22 inch Code 165

☐ Pair

spoke cover with special-motive

the price for special-spoke cover are net price and not discountable!

www.speichenschutz.com

motive no.: _____

please take the motive no. from the website www.speichenschutz.com

☐ 22" ☐ 20"

Code 2268

☐ Pair

www.rs-speichenschutz.com

motive no.: _____

please take the motive no. from the website www.rs-speichenschutz.com

☐ 22" ☐ 20"

Code 2267

☐ Pair

Example





Information about others spoke cover motives please take website

Technical changes and print mistakes are reserved.

DOAbility

Note:
the combination 22" and 6" at short frame (Code 2428 or 2429) is not possible !!!!

Castor

Castor with screw axle short, castor adapter M8 and Design fork black

☐ 4 inch Skater	solid	Code 530
☐ 4 inch grey	solid	Code 531
☐ 5 inch grey	solid	Code 532
☐ 6 inch grey	solid	Code 533
☐ 6 inch grey	pneumatic	Code 536

PU tyre in black

PU black

☐ 4 inch	Code 531S
☐ 5 inch	Code 532S
☐ 6 inch	Code 533S

Standard

Castor solid
with **plastic rim black**
and **black tyre**
(slim tyre)

☐ 4 inch	Code 534
☐ 5 inch	Code 535

☐ Pair

Design castor slim
black tyre (solid)

aluminium rim silver	**aluminium rim black**
☐ 4 inch Code 537	☐ 5 Zoll Code 538S
☐ 5 inch Code 538	

☐ Pair

Castor

Lightwheel Castor with **LED**

☐ 4 inch	Code 181LED
☐ 5 inch	Code 582LED

☐ Pair

Castor wheel more outside
and more forward
(forward 40mm / outside 38mm)

Code 2459



not possible at
Little AR with frame short

☐ Pair


Soft-aluminium-castor
black tyre (solid) **soft**

aluminium rim silver	**aluminium rim black**
☐ 4 inch Code 5394	☐ 4 inch Code 5394S
☐ 5 inch Code 5395	☐ 5 inch Code 5395S

☐ Pair

brake

Knee angle brake
offset with
short lever stroke



☐ Code 174K

Knee angle brake
"Pull to lock"

Pull to lock
the brake



☐ Code 274

Brake lever extension
straight, foldaway

(only for
CODE 174K)



Code 175

☐ Pair

Knee angle brake offset
with **locking pin** for locking






fixed in the **braked** position	fixed in the **braked** and **opened** position	fixed in the **opened** position
☐ Code 174Rast	☐ Code 174RastD	☐ Code 174RastO

Standard

brake

Knee angle brake 174K
switched through companion




Code 2478

Combination with 175 not possible
Combination with 2487 / 2489 not possible

☐ Pair

additonal Knee angle brake only for
companion, activated with
bowden cable

(Brake is mounted under the frame)



Code 2479

Combinaton with 175 not possible
Combination with 2487 / 2489 not possible

☐ Pair

accessories

Anchor plate system

anchor plates for transport in cars (ISO 7176-19)
Assembly kit for wheelchairs including mounting.
Little AR: from SD 36/38
Little SLT: possible

CRASH TESTED ISO 7176-19

Code 2490

☐ Set

note

Technical changes and print mistakes are reserved.

VD2362E-29

Little SLT — accessories

trakking extension **TRIKE** as steering and pushing aid

hook closure 

holder lateral at frame (for type d.23) with hook closure
casotr 10x2" black alu fork
arrangement with gear segment
with ramp as mounting aid

capacity up to **60kg** (wheelchair incl. payload)
size 1: frame straight SW 26cm - 28cm
frame abduction SW 24cm - 28cm
size 2: seat width 30cm - 34cm

CODE 9901
(ET-Nr.: T1010000)

☐ Piece

capacity up to **90kg** (wheelchair incl. payload)
size 3: seat width 36cm - 40cm
size 4: seat width 42cm - 46cm

CODE 9902

☐ Piece

Trike

trakking extension
as steering and pushing aid



Little AR

trakking extension **TRIKE** as steering and pushing aid

hook closure 

special holder lateral at frame (for type d.23) with hook closure
castor 10x2" black alu fork
Gestänge mit Zahnsegment
mit Rampe als Montagehilfe

capacity up to **60kg** (wheelchair incl. payload)
size 1: frame straight SW 26cm - 28cm
frame abduction SW 24cm - 28cm
size 2: seat width 30cm - 34cm

CODE 9906

☐ Piece

capacity up to **90kg** (wheelchair incl. payload)
size 3: seat width 36cm - 40cm
size 4: seat width 42cm - 46cm

CODE 9907

☐ Piece

TRIKE colour frame — Edition "....-Color"

Note:
the
TRIKE frame parts
are
made of aluminium

TRIKE- frame colour
powder coated

Colour : frame colour (see page 1)

- frame type right / left
- cross tube
- castor fork

Code 9908

without extra charge ☐

Edition "Duo-Color"

colour 1: _____
(castor fork)

colour 2: _____
(frame type right / left, cross tube)

Code 9909

☐


Castor fork with brake
for directional stability

Code 9909F

☐ Piece

Little SLT — accessories

trakking extension **TRIKE** as steering and pushing aid

quick release (clamp) 

holder quick release (clamp) at frame (for pipe Ø 23mm)
castor 10x2" black alu fork
arrangement of levers with gear segment
with ramp as mounting aid
with adjusting ring at frame (Ø 23mm)

capacity up to **60kg** (wheelchair incl. payload)
size 1: frame straight SW 26cm - 28cm
frame abduction SW 24cm - 28cm
size 2: seat width 30cm - 34cm

CODE 9921
(ET-Nr.: T1020000)

☐ Piece

capacity up to **90kg** (wheelchair incl. payload)
size 3: seat width 36cm - 40cm
size 4: seat width 42cm - 46cm

CODE 9922
(ET-Nr.: T1020020)

☐ Piece

Note:

Little AR: combination TRIKE with anchor plate system is **not** possible !!

Little SLT: combination TRIKE (hook) with anchor plate system is from frame long (from SD 34) possible !!
combination TRIKE (clamp) with anchor plate system is possible !!

Technical changes and print mistakes reserved.

DOAbility

Little

BEROLLKA bewegt

note

Technical changes and print mistakes are reserved.

DOAbility

General Terms and Conditions of Business (GTC)

Validity

All sales are subject to our delivery and payment conditions, without any exception. Orders and other agreements will not be valid unless confirmed by us in writing. We expressly object to the customer's purchase conditions.

Quotations and contract conclusion

Quotations are always subject to change, unless they are expressly limited in terms of time. Orders require a written confirmation. After confirmation, prices quoted in orders are fixed but we reserve the right to make adequate price adjustments in the event of price increases for which we are not responsible (e.g. tax increases, wage increases, price changes on the part of our suppliers). Cost estimates for customized products and repairs are always non-binding. We reserve the right to make technical modifications. The contract will not be concluded before the order is confirmed.

Shipment

To the extent the customer has not made any specification, shipment will be made in usual packaging using an adequate method. Insurance will only be taken out upon request and at the expense of the customer.
Delivery will be made carriage forward. We are free to select a method of shipment. If the customer has special requests, e.g. express delivery, express freight, etc., the additional cost will also be borne by the recipient. The goods will always be shipped at the buyer's risk. In case of damage during transport, the recipient must apply for an ascertainment of the facts with the respective forwarding agent.

Prices

All prices quoted by BEROLLKA will be quoted ex works, excluding any applicable statutory value added tax, public charges and levies such as potential customs tariffs, and insurance and packaging costs that will be charged separately.

Minimum quantity surcharge

If a small order is placed with a net value below 20 euros, an additional lump-sum minimum quantity surcharge of 15.00 euros will additionally be charged to cover administrative costs.

Delivery periods

The delivery dates and periods stated are subject to change. Any delay in delivery will not entitle the customer to withdraw from the contract. The assertion of damage claims is ruled out. BEROLLKA can make partial deliveries for justified reasons.

Wrong orders / return of merchandise

Any goods erroneously ordered will only be taken back provided we are responsible for such an error and the merchandise is in impeccable condition. In the event that an item erroneously ordered or merchandise that was returned is taken back in an individual case as a gesture of goodwill, attention must be paid to the following:

1. It is necessary to state the article number, the designation of the article, the date and number of the delivery note and the invoice number.
2. The return shipment must be carriage paid.
3. The merchandise must correspond to the article number(s) listed in the delivery note. In no event will any individual parts removed from an article be taken back.
4. The returned merchandise must be in impeccable condition.
5. Processing fees amounting to 10% of the net merchandise value not, however, less than €35.- will be charged for any instance where goods are taken back.

Deliveries effected more than 3 months earlier as well as sewing pieces and customized products cannot be exchanged. A copy of the delivery note or of the invoice will serve as proof in this regard.

Warranty and notices of defects

Excluding any further claims, the seller is liable for defects of the merchandise in that the seller will mend or replace all those parts free of charge that have become defective or unusable within 6 months due to provable defects of material or workmanship, and of which the seller has been notified in writing immediately after knowledge was gained of the defect. Should longer periods of time be required for statutory reasons, the warranty period is agreed to be 12 months.
The customer must grant the seller the required time and opportunity for effecting delivery of the replacement parts. Replaced parts will become the property of the seller. The buyer must check the merchandise without undue delay after it has been delivered by the seller and, if a defect is found, notify the seller of the defect without undue delay. Should the buyer fail to notify the seller, the merchandise will be deemed accepted unless a defect occurs that could not have been found when the merchandise was checked.
In any event in which a complaint is brought to the notice of BEROLLKA, BEROLLKA invariably has the right to inspect and check the item delivered that is the subject matter of said complaint, for which the customer will grant BEROLLKA the required time and opportunity. BEROLLKA can also demand that the customer sends the item delivered that is the subject matter of a complaint back to BEROLLKA at the expense of BEROLLKA. Should a notice of defects given by the customer turn out to be wilfully or gross negligently unfounded, the customer will be obliged to compensate BEROLLKA for any and all expenses incurred in this connection. BEROLLKA will be entitled to remedy warranty defects at its own discretion either by rectification of the defect or by replacement of the defective part or of the entire item delivered at no cost for the customer. Should the customer itself repair an item that gave rise to a complaint, or should the customer endeavour to rectify the defects itself without prior written approval of the seller, the seller will not grant any claims for reimbursement of expenses (such as transport, road, labour or material costs) or warranty claims.
The following are excluded from the liability for damages:
damage arising from natural wear and tear, damage due to negligent and/or improper treatment, or violent damage. A revocation of the purchase contract, a reduction of the purchase price or damage claims of any kind are ruled out. The warranty obligation will cease if repairs and/or modifications or replacements of individual parts are performed by other parties than the seller or the delivering company.

Reservation of title

With regard to resellers: The merchandise delivered will remain our property until any and all invoices – also for different deliveries – have been fully paid. Our merchandise that is hence still our property may only be resold within the ordinary course of business. Amounts receivable from any resale are assigned to us. We must be notified without undue delay and on the reseller's own initiative of any disposal of said merchandise outside of the ordinary course of business, e.g. through pledging our property or amounts receivable assigned to us.

Terms of payment

Any invoices are payable within 30 days from the date of the invoice strictly net cash. A 2% discount will be granted for payments received within 8 days provided any and all payment claims under earlier deliveries have been settled. We are entitled to charge default interest if the 30-day payment deadline is exceeded. The interest rate will be the customary market rate charged by major banks. We expressly reserve the right to effect delivery with cash on delivery.

Data Protection

BEROLLKA and the customer are responsible themselves for compliance with the respective data protection rules applicable to them. The customer must not provide BEROLLKA with personal data that is illegal or that the customer is not legally authorised to distribute.

The place of performance and the place of jurisdiction will be the place where the invoice is issued.

Heidelberg, Germany, will be the exclusive place of jurisdiction for any and all disputes arising under the contractual relationship. However, BEROLLKA is entitled to sue the customer at any other legal venue. These General Terms and Conditions of Business are governed by the laws of the Federal Republic of Germany, excluding of the United Nations Convention on Contracts for the International Sale of Goods of Goods (CISG).

Technical modifications and printing errors excepted. A current list of prices can be ordered at any time.